## ENGLISH

### DEVICE DESCRIPTION

The aspicio™ Toric Soft Hydrophobic Intraocular Lenses (TR60C, TR60Y) are foldable one-piece posterior chamber, Clear (UV absorbing) or Yellow (UV+blue light absorbing) optical implant lenses used for the replacement of the human crystalline lens in the visual correction of aphakia and pre-existing corneal astigmatism in adult patients, with or without presbyopia, who desire improved uncorrected distance vision, reduction of residual refractive cylinder, and increased spectacle independence for distance vision. The yellow aspicio™ Toric IOL contains ICARES Medicus' proprietary blue light filtering chromophore that filters light in a manner that approximates a young human crystalline lens in the 400-475 nm blue light wavelength range.

The aspicio™ Toric Soft Hydrophobic Intraocular Lens (TR60C, TR60Y) is a foldable posterior chamber, one-piece acrylic lens with square edge, a 6.0 mm optic, and an overall length of 13.0 mm. All models are manufactured at the following diopter power range: +4.0 D to +34.0 D in 0.5 diopter increments with an added cylindrical power of 1.0 to 6.0 D in 0.5 D increments as indicated by the axis orientation marks on the anterior optic surface to achieve visual correction of aphakia and pre-existing corneal astigmatism. Refractive index of the lens is 1.5.

The astigmatism to be corrected by the aspicio™ Toric IOL model TR60Y and TR60C shall be determined by keratometry and biometry data keeping in mind that the size and location of the surgical incision can affect the amount and axis of corneal astigmatism. To optimize the selection of the proper aspicio™ Toric IOL lens power, a web-based calculator ([www.icaresmedicus.com](http://www.icaresmedicus.com)) is provided for the surgeon. When the pre-operative keratometry and biometry data, incision location, and the surgeon's estimated surgically induced corneal astigmatism are entered, the web-based calculator will determine the appropriate TR60Y or TR60C Toric IOL lens spherical dioptric power, cylinder power, and axis of placement in the eye to be implanted.

### MODE OF ACTION

The aspicio™ Toric posterior chamber intraocular lenses are intended to be positioned in the posterior chamber of the eye, replacing the natural crystalline lens. This position allows the lens to function as a refractive medium in the correction of aphakia. The aspicio™ Toric IOLs have a toric aspheric optic with cylinder axis marks to denote the flat meridian (plus cylinder axis). Alignment of the toric IOL cylinder axis marks with the post-operative steep corneal meridian allows the lens to correct astigmatism. The toric aspheric optic reduces spherical aberration as compared to a standard spherical toric optic in an average eye.

### INDICATIONS

The aspicio™ Toric IOLs are indicated for primary implantation for the visual correction of pre-existing corneal astigmatism and aphakia in adult patients, with and without presbyopia, who desire improved uncorrected distance vision, reduction of residual refractive cylinder, and increased spectacle independence for distance vision following removal of a cataractous lens. The device is intended to be placed in the capsular bag by an ophthalmologist following successful circular tear anterior capsulotomy with a verified absence of radial tears.

### WARNINGS

*Physicians considering lens implantation should weigh the potential risk/benefit ratio for any circumstances that could increase complications or impact patient outcomes. This lens should not be implanted under the following conditions:*

- 1. If the posterior capsule is ruptured, if the zonules are damaged, or if a primary posterior capsulotomy is planned.*
- 2. The Tyvek pouch is found to be damaged or opened.*
- 3. Suspected microbial infection.*
- 4. Recurrent severe anterior or posterior segment inflammation or uveitis.*
- 5. Patients in whom the intraocular lens may affect the ability to observe, diagnose, or treat posterior segment disease.*
- 6. Surgical difficulties at the time of cataract extraction that might increase the potential for complications (e.g. persistent bleeding, significant iris damage, uncontrolled positive pressure, or significant vitreous prolapse or loss).*
- 7. A distorted eye due to previous trauma or developmental defect in which appropriate support of the IOL is not possible.*
- 8. Circumstances that would result in damage to the endothelium during implantation.*
- 9. Children under the age of 2 years are not suitable for intraocular lenses.*

### PRECAUTIONS

- 1. Do not resterilize the lens by any method.*
- 2. Do not store the lens at a temperature lower than 5°C or greater than 25°C.*
- 3. Do not reuse the lens. The lens is for single use only.*
- 4. Use only sterile intraocular irrigating solutions (such as BSS® or BSS PLUS®) to rinse/or soak lenses.*
- 5. Accurate keratometry and biometry in addition to the use of the Toric Calculator ([www.icaresmedicus.com](http://www.icaresmedicus.com)) are recommended to achieve optimal visual outcome.*
- 6. Handle lenses carefully to avoid damage to lens surface or haptics.*
- 7. Do not attempt to reshape haptics in any way.*
- 8. A high level of surgical skill is required for intraocular lens implantation. The surgeon should have observed and/or assisted in numerous implantations and successfully completed one or more courses on intraocular lens implantation before attempting to implant intraocular lens.*
- 9. Rotation of aspicio™ Toric IOLs away from their intended axis can reduce their astigmatic correction. Misalignment greater than 30° may increase postoperative refractive cylinder. If necessary, the lens repositioning should occur as early as possible prior to lens encapsulation. Some clinical cases suggest encapsulation is complete within four weeks of implantation.*

*10. Carefully remove all viscoelastic from both the anterior and posterior sides of the lens. Residual viscoelastic may allow the lens to rotate causing misalignment of the aspicio™ Toric IOL with the intended axis of placement.*

### CONTRAINDICATIONS

Patients with any of the following conditions may not be suitable candidates for an intraocular lens because the lens may exacerbate an existing condition, may interfere with diagnosis or treatment of a condition, or may pose an unreasonable risk to the patient's eyesight. Careful preoperative evaluation and sound clinical judgment should be used by the surgeon to decide the benefit/risk ration before implanting a lens in a patient with one or more of the following conditions.

Before Surgery:

1. Choroidal hemorrhage
2. Concomitant severe eye disease
3. Extremely shallow anterior chamber
4. Microphthalmos
5. Non-age-related cataract
6. Proliferative diabetic retinopathy (severe)
7. Severe corneal dystrophy
8. Severe optic nerve atrophy
9. Irregular corneal astigmatism
10. Medically uncontrolled glaucoma
11. Chronic severe uveitis
12. Diabetic retinopathy
13. Clinically significant macular/RPE changes

During Surgery:

1. Excessive vitreous loss
2. Capsulotomy by any technique other than a circular tear
3. The presence of radial tears known or suspected at the time of surgery
4. Situation in which the integrity of the circular tear cannot be confirmed by direct visualization
5. Cataract extraction by techniques other than phacoemulsification or liquefaction
6. Situation where the need for a large capsulotomy can be anticipated (e.g., diabetics, retinal detachment in the fellow eye, peripheral retinal pathology, etc.)
7. Posterior capsular rupture (preventing fixation of IOL)
8. Zonular damage (Preventing fixation of IOL)
9. Uncontrollable positive pressure
10. Significant anterior chamber hyphema

### COMPLICATIONS

The following lists potential complications accompanying cataract or implant surgery. Complications may include, but are not limited to the following:

Post-Operative Adverse Events:

1. Lens rotation
2. Visual distortion
3. Blurred vision at near distance

Cumulative Adverse Events:

1. Corneal endothelial damage
2. Infection (endophthalmitis)
3. Hyphema
4. Hypopyon
5. Lens Dislocation
6. Cystoid macular edema
7. Corneal edema
8. Pupillary block
9. Cyclitic membrane
10. Iris prolapse
11. Retinal detachment
12. Vitritis
13. Transient or persistent glaucoma
14. Secondary surgical intervention (excluding retinal detachment and posterior capsulotomy), include, but not limited to the following:
  - a. Iridectomy for pupillary block
  - b. Vitreous aspiration for pupillary block
  - c. Repositioning of lens
  - d. IOL removal for inflammation
  - e. IOL replacement
  - f. Wound leak repair

Persistent Adverse Events:

1. Corneal Stroma Edema
2. Cystoid Macular Edema
3. Iritis
4. Raised IOP requiring treatment

### SELECTION AND PLACEMENT OF THE aspicio™ TORIC IOL

The astigmatism to be corrected should be determined from keratometry and biometry data rather than refractive data since the presence of lenticular astigmatism in the crystalline lens to be removed may influence results. The size and location of the surgical incision may affect the amount and axis of corneal astigmatism. In order to optimize IOL selection and axis placement, ICARES provides a web-based tool ([www.icaresmedicus.com](http://www.icaresmedicus.com)) for the surgeon. Pre-operative keratometry and biometry data, incision location, and the surgeon's estimated surgically induced corneal astigmatism are used to determine the appropriate aspicio™ Toric IOL model, spherical equivalent lens power, and axis of placement in the eye.

For optimal results, the surgeon must ensure the correct placement and orientation of the lens within the capsular bag. The anterior surface of the IOL is marked with indentations (three at each end) at the haptic/optic junction that identify the flat meridian of the aspicio™ Toric IOL optic. These indentations form an imaginary line representing the plus cylinder axis (note: IOL cylinder steep meridian is 90° away). The aspicio™ Toric IOL cylinder axis marks should be aligned with the post-incision steep corneal meridian (intended axis of placement).

Prior to surgery the operative eye should be marked in the following manner:

With the patient sitting upright, precisely mark the twelve o'clock and/or the six o'clock position with a T marker, a surgical skin marker, or a marking pencil indicated for ophthalmic use. Using these marks as reference points, an axis marker can be used immediately prior to or during surgery to mark the axis of lens placement following the use of the web-based [www.icaresmedicus.com](http://www.icaresmedicus.com) to determine the optimal axis of placement.

After the lens is inserted, precisely align the axis marking indentations on the aspicio™ Toric IOL with the marked axis of lens placement. Carefully remove all viscoelastic from both the anterior and posterior sides of the lens. This may be accomplished by manipulating the IOL optic with the I/A tip and using standard irrigation/aspiration techniques to remove all viscoelastic from the eye. Bimanual techniques may be used, if preferred, to ensure removal of viscoelastic from behind the lens implant. Special care should be taken to ensure proper positioning of the aspicio™ Toric IOL at the intended axis following viscoelastic removal. Residual viscoelastic may allow the lens to rotate causing misalignment of the aspicio™ Toric IOL with the intended axis of placement.

Misalignment of the axis of the lens with the intended axis of placement may compromise its astigmatic correction. Such misalignment can result from inaccurate keratometry or marking of the cornea, inaccurate placement of the aspicio™ Toric IOL axis during surgery, an unanticipated surgically induced change in the cornea, or physical rotation of the aspicio™ Toric IOL after implantation. In order to minimize this effect, the surgeon should be careful to ensure that preoperative keratometry and biometry is accurate and that the IOL is properly oriented prior to the end of surgery.

### DIRECTIONS FOR USE

1. Examine the label on the lens box for proper lens model, diopter power and expiration date.
2. Open the lens box to remove the pouched lens and verify the lens case information (e.g. power, model and serial numbers) is consistent with the information on the outer box.
3. This device is sterile until the inner pouch is opened. Inspect the pouch carefully for tears, cuts, punctures or other signs that the pouch has been opened or damaged. DO NOT implant the IOL if the sterility has been compromised.
4. To remove the lens, open the undamaged pouch and transfer the case to a sterile environment. Carefully open the case to expose the lens.
5. To minimize the occurrence of marks on the lens due to handling, all instrumentation should be scrupulously clean. Any forceps used for lens handling must have round edges and smooth surfaces.
6. When removing the lens from the case, DO NOT grasp the optic area with forceps. The lens should be handled by the haptic portion only. Handle the lenses carefully to avoid damage to lens surfaces or haptics. DO NOT attempt to reshape haptics in any way.
7. Rinse the lens thoroughly using sterile intraocular irrigating solution such as BSS® or BSS PLUS® solution. DO NOT rinse the IOL in solutions other than sterile intraocular irrigating solution. Prior to insertion, the IOL should be carefully examined to ensure that particles have not adhered during handling.
8. ICARES Medicus, Inc. recommends using the lioli™ IOL Delivery System, pioli™ IOL Delivery System, or other equivalent commercial available delivery system. For lioli™ IOL Delivery System, LIOLI-18 only validated up to +20D hydrophobic intraocular lenses. LIOLI-22 and LIOLI-24 are validated up to +30D hydrophobic intraocular lenses. For pioli™ IOL Delivery System, PIOLI-D is validated up to +28D hydrophobic intraocular lenses. Failure to follow the instructions for use may result in patient injury. Sodium hyaluronate-based viscoelastic solution is recommended to be used with lioli™ and pioli™ IOL delivery systems. For the methods of loading the lens onto the injector, please refer to respective IFU of the IOL delivery system used during the surgical procedure.
9. There are various surgical procedures which can be utilized, and the surgeon should select a procedure which is appropriate for the patient. Surgeons should verify that appropriate instrumentation is available prior to surgery.
10. DO NOT reuse this IOL. The device is for single use only. Reuse could result in device contamination or damage to IOL or cartridge.
11. DO NOT RESTERILIZE.

### PATIENT REGISTRATION AND REPORTING

Contact local distributors regarding any reports of adverse events.

### HOW SUPPLIED

The aspicio™ Toric Soft Hydrophobic Intraocular Lens (TR60C, TR60Y) is supplied dry, in a lens tray packaged in a Tyvek peel pouch and terminally sterilized by ethylene oxide. The Lens must be opened only under aseptic conditions (See DIRECTIONS FOR USE above).

### EXPIRATION DATE

The packaged aspicio™ Toric Soft Hydrophobic Acrylic Intraocular Lens is sterile unless the Tyvek peel pouch is damaged or opened. The sterility expiration date is clearly indicated on the outside of lens package. DO NOT IMPLANT the aspicio™ Toric IOL after the expiration date. Any lens held after the expiration date should be returned to the distributor.

### TRANSMITTANCE COMPARISON

Model	Characteristics	UV Cutoff at 10% T
Yellow (TR60 Y)	UV-Absorber + blue light filter	≥ 384 nm
Non–Yellow (TR60 C)	UV-Absorber	≥ 382 nm

# aspicio™

## Soft Hydrophobic Intraocular Lens

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Rev. Z1

**ICARES™**  
MEDICUS

ICARES Medicus, Inc.  
4F, No.16, Sec.2, Sheng Yi Road, Zhubei, Hsinchu 30261, Taiwan, R.O.C.

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